





LABOMED, INC.

www.labomed.com
spectro@labomed.com



LB-1000 Comparison Microscope

LB-1000 Comparison Microscope is basically used for forensic science and police departments. It is also provided to these departments such as police training schools, banks, tax, print, coinage, archaeology and etc.

3.1×-120×Optical Magnification Range

12V 50W Halogen, adjustable, air-blast strong light source

12V 50W, adjustable, Φ10mm cold light source

This microscope illuminates by halogen lamp and focuses manually. The objectives are changed, one by one.