





# LABOMED, INC.

[www.labomed.com](http://www.labomed.com)  
[spectro@labomed.com](mailto:spectro@labomed.com)



## LB-1010 Digital Comparison Microscope

LB-1010 Digital Comparison Microscope has a built-in digital image analysis system, it is basically used for forensic science and police departments. It is also provided to these departments such as police training schools, banks, tax, print, coinage, archaeology and etc.

3Megapixel, supports USB/VGA/AV/DAI output, has camera connector

2.5×-120×Optical Magnification Range

This microscope comes with camera and software and illuminates by R, G, B LED light, with 50,000 hour lifetime, and changes color. It uses cold light and focuses manually. The objectives are changed, one by one.