






# LABOMED, INC.

[www.labomed.com](http://www.labomed.com)  
[spectro@labomed.com](mailto:spectro@labomed.com)



## LB-102 Digital Monocular Microscope

Magnification standard  $20\times \sim 135\times$ , the maximum can be  $1080\times$ , zoom magnification.

External light source (LED light) needed.

Measurement software: point-to-point straight-line distance, three-point Circle, radius, perimeter, area, angle, perimeter of the polygon, area measurement.

Alternative objective lens and eyepiece, so as to obtain a multiple magnification.

With modular structure, compact and complete accessories,

LB-102 is widely used in microelectronics manufacturing and testing industry, especially for high-precision test and inspection industry.