





# LABOMED, INC.

[www.labomed.com](http://www.labomed.com)  
[spectro@labomed.com](mailto:spectro@labomed.com)



## LB-120 Polarizing Monocular Microscope

0°-90° Rotatable Analyzer and when out of optical path, Single Polarizing Observation is Available. With Bertrand, Interferogram of Specimen can be observed clearly.  $\lambda$  Slip,  $1/4\lambda$  Slip and Quartz Wedge are Equipped  
Quadruple Nosepiece, Center Adjustable  
Polarizer can be Swung out of Optical Path  
6V/ 20W Halogen Lamp Illumination, Brightness Adjustable