



[www.labomed.com](http://www.labomed.com)  
[spectro@labomed.com](mailto:spectro@labomed.com)

## LB-122 Portable Monocular Inverted Biological Microscope for Educational Purposes

### Introduction

**LB-122 Portable Monocular Inverted Biological Microscope for Educational Purposes** is designed to observe liquid specimens and other filled samples. It is portable, easy to operate and has very good performance.

### Applications

**LB-122 Portable Inverted Biological Microscope for Educational Purposes** is ideally suited for school biological education. We can use it to observe all kinds of slides. We can also use it as a toy or other applications to examine coins, stamps, PC boards, insects, and other objects especially at the lower powers.

### Technical Specifications

Eyeiece:	Eyeiece 10×/18
Achromatic Objective:	4× NA0.08/WD16 10× NA0.25/WD3.6 40× NA0.60/WD0.6
Nosepiece:	Triple Nosepiece
Stage:	Plain Stage with Slide Clips
Condenser & Illumination:	Bright field, LED Illumination(AA Battery X2)
Moving Range:	2.6mm

