





# LABOMED, INC.

[www.labomed.com](http://www.labomed.com)  
[spectro@labomed.com](mailto:spectro@labomed.com)



## LB-1261 Compound LCD Digital Biological Microscope with Infinity Optical System (5.0MP)

Built-in 5 million-pixel digital camera, image data can be easily retained without external computers and other equipment, they can improve the efficiency of research and analysis.

8-inch digital LCD screen, high definition and bright colors, easy for people to share. Large storage capacity, with SD card, print pictures instantly without being connected to PC.