




LABOMED, INC.

www.labomed.com
spectro@labomed.com



LB-1262 Compound LCD Digital Biological Microscope with Infinite Optical System, LED Illumination and Software

The LB-1262 LCD digital microscope is a new generation microscope that with an embedded operation system WinCE 5.0, it can link mouse or keyboard, just like a micro computer.

It supports 100M Ethernet network, WI-FI wireless net and can communicate with remote PC.

With 8.4 Inch TFT touching screen, friendly man-machine interface, multiple peripheral interfaces and powerful processing function, make it very easy to operate.

It can be used in various optical micro fields, such as teaching, biological and medical researching, electronic components inspection and so on.

The screen can be separated from the microscope and placed anywhere you like.