




LABOMED, INC.

www.labomed.com
spectro@labomed.com



LB-1601 Digital LCD Metallurgical Measuring Microscope with Infinity Optical System (5.0MP)

Special designed metallurgical analysis optical system and an 8" LCD screen, images and videos can be share with other people.

5 Megapixels with real time imaging display, images and video recording and re-playing function.

X/Y digital micrometer and moving stage for accurate measuring.

12V DC power supply, safety and reliability of low voltage (Input Voltage 90-240V).

USB Cable for transferring data to a PC.

SD Card slot for storage, printing, optional SD cards (Maximal Storage Capacity is 4G).

Multi-function top LED ring illuminating system, four-region can be controlled separately.