




# LABOMED, INC.

[www.labomed.com](http://www.labomed.com)  
[spectro@labomed.com](mailto:spectro@labomed.com)



## LB-1910 Surgical Operation Microscope

All the lenses are multi-coated, mildewproof and anti-reflective.

Optional objectives: f250/f300/f350f/f400.

Desktop stand can be chosen to make the device more portable. It can also be customized according to various operation requirements.

Foot control focus, 3 gears magnification, sharp images and comfortable observation, can meet the needs of cataracts and other eye surgeries.

With apochromatic technology, focal points of different wavelengths light through the lens are very close, the operators' vision is more clearly.

OPTIONAL: Video Devices (1/2 or 1/4" CCD Camera, CCD Adapter, Beam Splitter and Connecting Cable)