




# LABOMED, INC.

[www.labomed.com](http://www.labomed.com)  
[spectro@labomed.com](mailto:spectro@labomed.com)



## LB-1940 Surgical Operation Microscope

ZOOM magnification system,  $4\times \sim 24\times$  zoom magnification, horizontal movement of foot control and fine focusing.

Compact structure and complete functions.

Germany imported Schott optics ensure the sharp image and even field of view.

Cold light source fiber illumination does no harm to tissues.

OPTIONAL: Video Devices (1/2 or 1/4" CCD Camera, CCD Adapter, Beam Splitter and Connecting Cable)