




LABOMED, INC.

www.labomed.com
spectro@labomed.com



LB-352 Stereo Zoom Digital Trinocular Microscope With Greenough Optical System and 5MP Color Digital Camera

Updated optical design, large depth of field, stereoscopic effect.

Ergonomic design, sharp image display, wide viewing field and easy to operate.

Providing a wide viewing field of $\Phi 23\text{mm}$ with the zoom ratio of 1:7.6.

Transmitting & Reflecting LED Illumination.

Ideal instrument in academic, medical and industrial fields.