



www.labomed.com
spectro@labomed.com

LB-933 Multi-Head Microscope (for 5 Users) (Anti-Mold)

Introduction

LB-933 Multi-Head Microscopes (for 5 Users) are equipped with multi-heads for more persons observing at the same time. The optical system is high quality and reliable.

Features

Sharp image display with excellent optical system.

Small space occupancy with integral stand design, low environment requirement with anti-mold technology.

Users friendly and comfortable operation with ergonomic design.

Applications

LB-933 Multi-Head Microscopes (for 5 Users) are mainly used in medical teaching and biological teaching. They also can be used for biological analysis for more experts at the same time. They can be widely used in medical and sanitary establishments, laboratories, institutes, colleges and universities.

Technical Specifications

| | | |
|---------------|---|-------|
| Viewing Head: | Sliding Binocular Head, inclined at 45° | 5Pcs |
| Eyepiece: | Wide Field Eyepiece WF10× | 10Pcs |
| Nosepiece: | Quadruple Nosepiece | |
| Objective: | Achromatic Objective 4×, 10×, 40×, 100× | |
| Condenser: | Abbe Condenser NA1.25 | |
| Focusing: | Coaxial Coarse & Fine Focus Adjustment, Fine Division 0.002mm | |
| Stage: | Plain Stage 140×140mm/ 75×50mm | |
| Illumination: | Halogen Lamp 6V/ 20W Halogen Lamp 24V/ 100W | |
| Pointer: | Green LED Pointer, Brightness Adjustable | |

