




LABOMED, INC.

www.labomed.com
spectro@labomed.com



LB-228 Binocular Digital Microscope with Compensation Free Binocular Head, Six Plan Achromatic Objectives ($4\times$, $10\times$, $20\times$, $40\times$, $60\times$, $100\times$) Camera and Software (1.3MP) and LED Illumination

Integrated digital microscope, with 1.3 Megapixel CMOS image sensor, the output images have high Signal to Noise Ratio;

High definition images and videos, high frame rate;

Professional software of real-time review, take photos and videos, measurement. Easy to operate;

LED light illumination, saves energy;

Compact and flexible, ideally suited for desktop, laboratory worktable;