




LABOMED, INC.

www.labomed.com
spectro@labomed.com



LB-286 Motorized Auto-Focus Digital Microscope

With four motors drive, objective nosepiece, stage with X, Y, Z directions can be controlled manually or motorized. With advanced digital illumination adjustment technology, brightness of halogen lamp has been divided into 100 steps.

3.2 Megapixels camera system, USB 2.0 plugs and play; real-time preview is up to 30 frames/ sec.

High precision stage, quickly capture the clearest images by auto-focus. Automatically slide scanning, recognizing are controlled by the computer.

With built-in dynamic focusing system, realize the laboratory microscopy automation. With network remote control system, it is widely used in medical diagnosis, pathological analysis and real-time teaching etc.